

Inventors: Lockridge and Watkins  
Serial No.: 09/748,739  
Filed: December 26, 2000  
Page 4

ay  
consists of 13 base-pair repeats separated by an 8 base-pair  
spacer: 5'-GAAGTTCCTATTC[TCTAGAAA]GTATAGGAACTTC-3' (SEQ ID NO:  
24). Briefly, variant libraries corresponding to the region of  
butyrylcholinesterase corresponding to amino acids 277-289 (SEQ  
ID NO: 13) of butyrylcholinesterase (shown as region 5 in Table  
2) were transfected into mammalian cells using flp recombinase  
and the 293T cell line. Table 5 shows the butyrylcholinesterase  
variants S285G, P285Q and P285S that were identified and  
characterized using the methods described herein utilizing Flp  
recombinase and the 293T human cell line.

---

REMARKS

The specification has been amended in four places to  
insert sequence identifiers. Concurrently with this amendment,  
Applicants are filing a paper copy and computer readable form of  
the sequence listing along with the appropriate statement under  
Rule 821 subsections (g) and (f).

Applicants have set forth above the amendment to the  
specification in clean form as required under 37 C.F.R. §  
1.121(b)(1)(i) and (ii). Applicants also attach Appendix A with  
marked up amendments indicated with brackets and underlining as  
required under 37 C.F.R. § 1.121(b)(1)(iii). Entry of the above-  
proposed amendments to the specification is respectfully  
requested.


Inventors: Lockridge and Watkins  
Serial No.: 09/748,739  
Filed: December 26, 2000  
Page 5

**CONCLUSION**

In light of the Amendments and Remarks herein,  
Applicants submit that the claims are now in condition for  
allowance and respectfully request a notice to this effect.  
Should the Examiner have any questions, he is invited to call  
Cathryn Campbell or the undersigned attorney.

Respectfully submitted,

May 29, 2002  
Date

  
Astrid R. Spain  
Registration No. 47,956  
Telephone: (858) 535-9001  
Facsimile: (858) 535-8949

CAMPBELL & FLORES LLP  
4370 La Jolla Village Drive  
7<sup>th</sup> Floor  
San Diego, California 92122  
**USPTO CUSTOMER NO. 23601**

Attachements: Appendix A - marked up version of amendments

Inventors: Lockridge and Watkins  
Serial No.: 09/748,739  
Filed: December 26, 2000

#### APPENDIX A

At page 25, Table 1:

Human BchE ( <u>SEQ ID NO: 25</u> )	act cct ttg tca gta
S287G ( <u>SEQ ID NO: 26</u> )	act cct ttg <b>ggt</b> gta
P285Q ( <u>SEQ ID NO: 27</u> )	act <b>cag</b> ttg tca gta
P285S ( <u>SEQ ID NO: 28</u> )	act <b>tcg</b> ttg tca gta

The paragraph spanning pages 43 and 44:

One approach for targeting variant or heterologous nucleic acids to a single site in the genome uses Cre recombinase to target insertion of exogenous DNA into the eukaryotic genome at a site containing a site specific recombination sequence (Sauer and Henderson, Proc. Natl. Acad. Sci. USA, 85:5166-5170 (1988); Fukushima and Sauer, Proc. Natl. Acad. Sci. U.S.A. 89:7905-7909 (1992); Bethke and Sauer, Nuc. Acids Res., 25:2828-2834 (1997)). In addition to Cre recombinase, Flp recombinase can also be used to target insertion of exogenous DNA into a particular site in the genome (Dymecki, Proc. Natl. Acad. Sci. U.S.A. 93:6191-6196 (1996)). The target site for Flp recombinase consists of 13 base-pair repeats separated by an 8 base-pair spacer: 5'-GAAGTTCCTATTC[TCTAGAAA]GTATAGGAACTTC-3' (SEQ ID NO: 24). As described herein, the butyrylcholinesterases designated SEQ ID NOS: 4, 6, and 8, were obtained by transfection of variant libraries corresponding to region 5 of human butyrylcholinesterase (see, Table 2) into mammalian cells using Flp recombinase and the human 293T cell line. It is understood

Inventors: Lockridge and Watkins  
Serial No.: 09/748,739  
Filed: December 26, 2000

that any combination of site-specific recombinase and corresponding recombination site can be used in methods of the invention to target a nucleic acid to a particular site in the genome.

At pages 58 and 59, please amend the paragraph spanning pages 58 and 59 to read as follows:

Thirty-four variants were prepared using PCR-site directed mutagenesis of human butyrylcholinesterase DNA performed utilizing Pfu polymerase (Stratagene, La Jolla, CA). Three oligonucleotide primers were used to perform the mutagenesis. The mutagenesis primers were used at the same time as a general primer such as the SP6 promoter sequencing primer (MBI Fermentas, Amherst, NY) to amplify one end of the butyrylcholinesterase cDNA. The following primers were used to prepare the A328W mutant: A328W antisense 5' ATAGACTAAAAACCATGTCCCTTCATC 3' (SEQ ID NO: 29); T7 old sense 5' TAATACGACTCACTATAGGG 3' (SEQ ID NO: 30); and SP6 antisense 5' ATTTAGGTGACACTATAG 3' (SEQ ID NO: 31). The A328W primer spans 27 nucleotides and contains the A328W mutation in the middle of the primer. The PCR reaction products (megaprimers) were cleaned on QuiaQuick PCR (Qiagen, Santa Clarita, CA) according to the manufacturer's protocol to remove excess primers. The cleaned megaprimers were extended in a second PCR reaction to generate the complete 1.8 kb coding sequence of each of the 34 variants.

Inventors: Lockridge and Watkins  
Serial No.: 09/748,739  
Filed: December 26, 2000

At page 77, first paragraph:

As shown in Table 5, several cell lines as well as other transfection methods were also characterized. As disclosed herein, Flp recombinase also can be used to target insertion of exogenous DNA into a particular site in the genome as described by Dymecki, supra, 1996. The target site for Flp recombinase consists of 13 base-pair repeats separated by an 8 base-pair spacer: 5'-GAAGTTCCTATTC[TCTAGAAA]GTATAGGAAGCTTC-3' (SEQ ID NO: 24). Briefly, variant libraries corresponding to the region of butyrylcholinesterase corresponding to amino acids 277-289 (SEQ ID NO: 13) of butyrylcholinesterase (shown as region 5 in Table 2) were transfected into mammalian cells using flp recombinase and the 293T cell line. Table 5 shows the butyrylcholinesterase variants S285G, P285Q and P285S that were identified and characterized using the methods described herein utilizing Flp recombinase and the 293T human cell line.